


The Grid-500 is a grid structure for relocating events on a glass coverslip. It provides 400 distinguishable observation squares of 500  $\mu\text{m}$  edge length. The grid is clearly visible by microscopy and imprinted into a microscopy coverslip. The outer dimensions and other parameters are identical to ibidi  $\mu$ -Dishes.

### Geometry of the Grid-500



- A to K (J not used) and 1 to 10
- A to K (J not used) and 11 to 20
- L to U and 1 to 10
- L to U and 11 to 20

#### Geometry of the Grid-500

Number of squares	400
Repeat distance	500 $\mu\text{m}$
Groove width	20 $\mu\text{m}$ ( $\pm 5 \mu\text{m}$ )
Groove depth	< 5 $\mu\text{m}$



4 x 10 x 10 squares




Microscopic images of the grid. Left: 4x objective lens brightfield without cells. Right: 10x objective lens phase contrast with rat fibroblast cells.

The four major squares are separated in 10 x 10 observation fields and indicated by letters and numbers ranging from:


### Characteristics of the Grid

The Grid-500 is made of small grooves that are imprinted into a microscopy coverslip. The structure is imprinted on the side on which cells are growing. Cells and grid are in one focal plane. There is no reported effect on cell growth, coating protocols, or surface properties. Proliferation and cell behavior is comparable to standard non-gridded glass coverslips. Washing steps (e.g. with PBS) before cell seeding can remove glass dust which is advantageous for direct cell growth on the surface. Please also refer to the instructions of ibidi  $\mu$ -Dishes glass bottom for information on surfaces, coatings, and cell seeding.


The grooves are 20  $\mu\text{m}$  ( $\pm 5 \mu\text{m}$ ) wide and approximately 5  $\mu\text{m}$  deep. Cells can grow in the grooves as well. We recommend using objective lenses up to 20x. Anyhow, the optical quality meets the requirements of 63x and 100x oil objective lenses as well (glass coverslip bottom, No. 1.5).




**μ-Dish <sup>35 mm, high</sup> Grid family**

 μ-Dish <sup>35 mm, high</sup> with Grid-500



Ordering number	Treatment or Coating	Characteristics
81166	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
81161	uncoated, sterile	hydrophobic

 μ-Dish <sup>35 mm, low</sup> with Grid-500


Ordering number	Treatment or Coating	Characteristics
80156	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
80151	uncoated, sterile	hydrophobic

 μ-Dish <sup>35 mm, high</sup> glass bottom with Grid-500


Ordering number	Treatment or Coating	Characteristics
81168	glass bottom, sterile	uncoated glass coverslip

 μ-Dish <sup>35 mm, high</sup> glass bottom with Grid-50


Ordering number	Treatment or Coating	Characteristics
81148	glass bottom, sterile	uncoated glass coverslip

**For research use only!**

Further technical specifications can be found at [www.ibidi.com](http://www.ibidi.com). For questions and suggestions please contact us by mail [info@ibidi.de](mailto:info@ibidi.de) or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany. ©ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.